

Position Pro II Insert

- Corrects “hammocking”
- Reduces pressure
- Discourages sacral sitting & sliding
- Increases sitting stability
- Discourages hip internal rotation and adduction
- Increases sitting tolerance

Matrix & Position Pro II

1. Unzip cover. Locate pocket on the **INSIDE** and **BOTTOM** of the cushion cover (the pocket is comprised of a stiff vinyl fabric).
2. Slide *Position Pro Insert* into pocket with the **label side up** and **arrows pointing toward the front** of the cushion. Slide *Position Pro* forward **until it rests against the front** of the pocket.
3. Close the pocket by **pressing together the hook & loop** (Velcro).


NOTE: *if your cover does not have hook & loop (Velcro) across the pocket opening, use the enclosed adhesive backed hook and loop to secure the Position Pro insert so that it is resting all the way forward in the pocket.*

4. Re-zip cushion cover closed.

Silver Care Basic & Position Pro II

1. Locate pocket on the **OUTSIDE** and **BOTTOM** of the cushion cover.
2. Slide *Position Pro Insert* into pocket with the **label side up** and **arrows pointing toward the front** of the cushion. Slide *Position Pro* forward **until it rests against the front** of the pocket.
3. Close the pocket by **pressing together the hook & loop** (Velcro).

NOTE: *if your cover does not have hook & loop (Velcro) across the pocket opening, use the enclosed adhesive backed hook and loop to secure the Position Pro insert so that it is resting all the way forward in the pocket.*



Remove *Position Pro* insert before washing the cushion.